

HOST

(Hierarchy Of Speech Targets)

Software Manual

1. Program Description

HOST (Hierarchy Of Speech Targets) is a universal open-source software which can be used remotely in telemedicine that will allow diagnosis of speech impaired children all around the world. Especially, in remote location where speech-language pathologists are not available, children born with Cleft Palate who are underserved in the rural area.

The desktop application was developed by the funding of National Institute of Dental and Craniofacial Research (Grant # 1UG3DE030869).

Principal Investigator:

Dr. Sue Ann Lee, Texas Tech University Health Sciences Center.

Subcontractors:

Dr. Jo Woon Chong, Texas Tech University.

Dr. Francisco Ortega, Texas Tech University.

Developers:

Kamrul H Foysal, Texas Tech University.

Bipasha Kundu, Texas Tech University.

Yifan Li, Texas Tech University.

2. Program Requirements

HOST (Hierarchy Of Speech Targets) software is developed in python 3 platform with tkinter, panda, numpy, openpyxl and other packages. The beta version (v1.0) requires Windows 7 to 11, with 2GB RAM (4GB preferable).

3. *Installation*


The software is sent in whole package as a zip file. Therefore, no installation is necessary. The user may use the software by only clicking the software (.exe file).

Alternatively, the software can be used with a shortcut which may be convenient to the user.


4. Starting the program:





1. Double-Click on the icon to start the program.




2. The Start-Screen will appear.


3. Button session type  shows the options for the sessions.



4. Button About  shows the description and developer information of the software.




5. Button Start Session  starts the program.

6. Button Exit Program  exits the program.

5. Software Flowchart:

The software is a complex web of logics meticulously crafted to evaluate the speech profile of a participant. The simplest flowchart of the software is shown below:




Pre-Treatment: Assessment of the participant before the traditional biofeedback evaluation

Traditional Biofeedback: Traditional biofeedback evaluation process incorporates sound level, syllable level, word level, phrase level and sentence level prompts for evaluation of the participants.

Post-Treatment: Assessment of the participant after the traditional biofeedback evaluation


6. Profile Creation:

User profile is created using the window below. After the program starts, this window will appear to enlist the participant.




The "Speech Therapy Application" window features a list box on the left labeled "Patient Name" with an empty space and a vertical scrollbar. To the right, under "Type of Session", are three radio button options: "Pre Treatment", "Traditional Biofeedback", and "Post Treatment". At the bottom left are "Add Client" and "Remove Client" buttons, and at the bottom center is a "Start" button.

6.1 Add Client:




The "Profile Creation" window contains the following fields: "Patient Name" (text box), "Year" (text box), "Month" (text box), "Day" (text box), a "Calculate" button, "Age" (text box), "Gender" with "Male" and "Female" radio buttons, a "Choose a Date" dropdown menu showing "5/25/22", a "Comments" text box, and "Ok" and "Cancel" buttons at the bottom.



The "Profile Creation" window is shown with the following filled values: "Patient Name" is "John Doe", "Year" is "2016", "Month" is "01", "Day" is "01", "Age" is "6", "Gender" has "Male" selected, "Choose a Date" is "5/25/22", and "Comments" is "Bilingual". The "Calculate" button is visible between the date and age fields.


Pressing the OK button includes the name of the participant in the list box.



Increasing number of participants will additionally show a scroll bar to access the participant.

6.2 Remove Client:

Selecting a participant and clicking the remove client button removes the entire profile of the participant.




6.3 Multiple Visit:

For multiple visits, the user do not need to create profiles for the same participants. Just clicking on the name of the participant will allow the user to conduct the study with further input.

6.4 Types of Session:

The radio buttons include the types of sessions.




6.5 Saving the Data:

When the user adds a participant, a folder is created with the first visit as visit 0. The program also creates 3 separate excel files automatically, each with individual info. 'Database_subjects.xlsx' saves the patients coming to the evaluation process. Under the participant's name, 'Subject_name_database.xlsx' is created with the subject information. The 'Subject_name_database.xlsx' saves the visit information.

7. Pre-Treatment:


7.1. **Flowchart:** The pre-Treatment flowchart (simplified user version) is depicted below.



7.2. Starting Pre-Treatment Session:

Selecting the Pre-Treatment Radio button starts the pretreatment session.

Clicking the **start button** with the pretreatment option will prompt user to select 4 different sounds




from the check box window below:

.

The **Show Result** button will start the evaluation process.

7.3. Counter:

At the right top of the window, the counter will show the number of prompts.



Pressing '1' for 'Correct' and '0' for 'Incorrect' will be used for evaluation process.

This process selects 16 words of each sound with visual cues and shows all 64 sounds randomly. At the end of the evaluation process, the screen shows the user the result on the screen.


Pressing other buttons other than 1s and 0s will not change the evaluation process until the user presses the correct button.

7.4. Saving the data:

The evaluation data is saved in the 'visit 0' folder under the participant's name in the 'Speech Therapy Application' program of the software mother folder. The excel file named 'Pre-Treatment Assessment.xlsx' saves the assessment information.


8. Traditional Biofeedback:

8.1. **Flowchart:** The Traditional BioFeedback flowchart (simplified user version) is depicted below.



8.2. Starting Traditional Biofeedback Session:

Selecting the Pre-Treatment Radio button starts the pretreatment session.



Clicking the **Start** button with the Traditional Biofeedback option will prompt user to the window below:

Target Speech Sounds

Visit summary

Visit No	Sound 1	Level sound 1	Sound 2	Level sound 2
----------	---------	---------------	---------	---------------

stop: ☐ b ☐ p ☐ t ☐ d ☐ k ☐ g

fricative: ☐ f ☐ v ☐ s ☐ z ☐ θ ☐ ç ☐ h

affricate: ☐ tʃ ☐ dʒ

nasal: ☐ m ☐ n

others: ☐ w ☐ j ☐ l ☐ r

Vowel Context: ☐ high ☐ Mid/ Low

Level Accuracy: ☐ 12 ☐ 24 ☐ 36 ☐ 48 ☐ 60

Level:

Show Result

The window includes the visit summary list box, Sounds, Vowel Context, Level of Accuracy and Level options.

For the first visit, the list box will show empty. Selecting a sound, selecting level of accuracy and putting level number and clicking the Show Result button starts the evaluation process.

Visit summary

Visit No	Sound 1	Level sound 1	Sound 2	Level sound 2
1	w	45	nan	nan

Initially, the user will put level 1. However, it is up to the user from which level s/he wants to start the evaluation process.

8.3. Visual Feedback:

8.3.1. Block and Level completion feedback:

An animation is shown after the participant completes a block or a level. This is to engage the participants.


8.3.2. Word level feedback:

Words are also visually shown with cues.

8.4. Feedbacks:

Feedback	Response
Correct	1
Incorrect	0
Good Job	2
Needs to improve	9

The feedbacks will be shown at the end of the session on the screen.




8.5. Counter:


At the right top of the window, the counter will show the number of prompts in the level.

8.6. Time Constraint:

The evaluation process will automatically shut down for that sound at 10 minutes. At 7 minutes.30 seconds, the program prompts a warning message.



After the window shuts down, the user can start again with the same sound or a different sound.




8.7. Saving the Data:

When the user selects a sound, the whole evaluation prompts for 150 levels are created and saved in result_X.xlsx file. For two individual sounds, two separate files (result1.xlsx and result_X.xlsx) will be created.

Each of these files will save the prompts, responses (1s and 0s) and feedbacks (2s and 9s).


9. Post-Treatment:

9.1. **Flowchart:** The Post-Treatment flowchart (simplified user version) is depicted below.




9.2. Starting Pre-Treatment Session:

Selecting the Post-Treatment Radio button starts the post-treatment session.



Clicking the start button with the post- treatment option will prompt user to select 4 different sounds

Clicking the **start button** with the pretreatment option will prompt user to select 4 different sounds



from the check box window below:

Target Speech Sounds

stop: ☐ b ☐ p ☐ t ☐ d ☐ k ☐ g

fricative: ☐ f ☐ v ☐ s ☐ z ☐ θ ☐ ç ☐ h

affricate: ☐ tʃ ☐ dʒ

nasal: ☐ m ☒ n


others: ☐ w ☒ j ☒ l ☒ r

Show Result

The **Show Result** button will start the evaluation process.

8.3. Counter:

At the right top of the window, the counter will show the number of prompts.



Pressing '1' for 'Correct' and '0' for 'Incorrect' will be used for evaluation process.

This process selects 16 words of each sound with visual cues and shows all 64 sounds randomly. At the end of the evaluation process, the screen shows the user the result on the screen.

Pressing other buttons other than **1s** and **0s** will not change the evaluation process until the user presses the correct button.

8.4. Saving the data:

The evaluation data is saved in the 'visit n' (n= last visit) folder under the participant's name in the 'Speech Therapy Application' program of the software mother folder.

Customizing the software:

There is option for customizing the software (time constrain, segment customization, feedback response and title/prompt). If any specific functionality is required, the developers can be contacted to customize the program to the users' choices.

Reference:

Jo Woon Chong, Ph.D.

Assistant Professor
Electrical and Computer Engineering
Email: J.Chong@ttu.edu
Phone: 806.834.8392
Fax: 806.742.1245

Sue Ann S. Lee, Ph.D., CCC-SLP

Professor
Dept of Speech-Language-Hearing Sciences
Email: sueann.lee@ttuhsc.edu
806-743-9051

Kwanghee Jung, Ph.D.

Assistant Professor
Educational Psychology, Leadership, & Counseling
Email: kwanghee.jung@ttu.edu
Phone: (806) 834-0023

Software Update log

May 31st, 2022

Kamrul H Foysal, Bipasha Kundu, Yifan Li.

Updated Speech Therapy Intervension for the version 1.0 (Beta).

June 8th, 2022

Kamrul H Foysal, Bipasha Kundu, Yifan Li.

Updated Speech Therapy Intervension for the version 2.0 (Beta).